



**Alcohol and Drug  
testing specialists**

# **SoToxa™** Oral Fluid Mobile Test System

**Rapid screening for drugs anywhere, anytime**



**Can test in a wide  
range of environments**



**Rapid screening for drugs of abuse in oral fluid in just 3 steps**

**THC (Dagga) sensitivity of 25ng/mL**

**Robust and able to test in a wide range of environments**

**Enhanced data management capabilities**

**Compact and portable**

**Ensuring a safer working environment for over 40 years**

**visit [www.alcosafe.co.za](http://www.alcosafe.co.za) for more information**

# SoToxa™ Oral Fluid Mobile Test System

We've designed and developed the So Toxa™ Oral Fluid Mobile Test System to meet the high performance demands of law enforcement agencies, drug treatment centres and workplace environments. The result is a drug testing solution which combines speed, ease of use and reliability.


## Key features

- Rapid oral fluid collection
- Increased sensitivity to THC (dagga)
- Easy to use
- Minimal chance for operator error
- Portable, compact and lightweight
- CE marked
- Full colour screen, readable day and night
- Stores 10,000 results
- Software to manage test data
- Instant printout


Drug Class	Cut-off (ng/ml)
Amphetamine	50
Benzodiazepines	20
Cocaine	30
Methamphetamine	50
Opiates	40
THC	25

## Drug testing as easy as


**A**  
Insert test cartridge into analyser



**B**  
Collect oral fluid sample



**C**  
Insert collection device into test cartridge





**The printouts have proven to be valuable evidence in disciplinary hearings**

## SoToxa™ Data Manager

- Download test data to a PC
- Compile reports to show trends, positivity rates and census information
- Add donor and tester names to results data
- Backup and store test data
- Configure settings and apply updates to the SoToxa™ Oral Fluid MobileTest System

